



Urinalysis reagent for visual or machine reading

AUTION Sticks

AUTION Sticks 10V



History of excellence ARKRAY's know-how and expertise has resulted in the development of AUTION Sticks, the most advanced urine test strips.

Urinalysis reagent for visual or machine reading

AUTION Sticks

AUTION Sticks 10V

Related products

Compact Urine Analyzer

Pocketchem™ UA

PU-4010

Code: 102678



Versatile Urine Analysis System

AUTION ELEVEN™

AE-4020

Code: 102676



Fully Automated Urine Analyzer

AUTION MAX™

AX-4030

Code: 102674



AUTION Sticks provide both visual and automatic measurement with ARKRAY'S exclusive automated analyzer.

Specification

	Measurement principle	Measurement range	Visual measurement time
Glu	Glucose oxidase reaction	50-1000 mg/dl	60 sec.
PRO	Protein-error reaction	15-1000 mg/dl	60 sec.
BIL	Azo-coupling reaction	0.5-6.0 mg/dl	60 sec.
URO	Azo-coupling reaction	2-8 mg/dl	60 sec.
pH	pH indicator	pH 5-9	60 sec.
S.G.	Cation extraction	S.G. 1.000-1.030	60 sec.
BLD	Activity measurement of pseudoperoxidase in hemoglobin	hemoglobin 0.03-1.0mg/dl	60 sec.
KET	Sodium nitroprusside method	5-150 mg/dl	60 sec.
NIT	Griess reaction	nitrite 0.08-0.5 mg/dl	60 sec.
LEU	Activity measurement of esterase in leukocytes	25-500 Leu/μl	90 sec.

ARKRAY GLOBAL BUSINESS, INC.

Yousuien-nai, 59 Gansuin-cho,
Kamigyo-ku, Kyoto 602-0008, JAPAN

TEL 81-75-662-8979 FAX +81-75-431-1202

