

Advanced Blood Glucose Test Meter

GLUCOCARD™ 01-mini plus

GT-1970



More advanced SMBG meter which successfully benefits the user with no coding, top specifications, and stylish design.

Advanced Blood Glucose Test Meter

GLUCOCARD™ 01-mini plus

GT-1970



Design

Impressive & Stylish design

Operation

Comprehensive display
New user friendly functions

Small Sample

0.3uL micro sampling

Maximum flexibility

Auto coding
Short fill detection
Strip discard button
Test strip port light
Backlight LCD
External output



Test strip for GLUCOCARD 01-mini plus

Specification

Test	Glucose in whole blood	Memory	500 test results
Sample Volume	Approx. 0.3uL	Operating Environment	Temperature range 10 to 40°C (50 to 104°F) Humidity 20 to 80% RH (Without condensation)
Test Strip	GLUCOCARD 01-SENSOR	Dimensions	102x33x18.5mm
Measuring Range	10 TO 600mg/dL (0.6 to 33.3 mmol/L)	Weight	Approx. 30g
Measuring Time	7 seconds		
Power source	2 X 3-volt Lithium battery (CR2032)		

※ Design and specifications may be changed without prior notice

ARKRAY GLOBAL BUSINESS, INC.

Kyoto Miyuki Bldg. 10F, 689 Takanna-cho,
Nakagyo-ku, Kyoto 604-8153, JAPAN

TEL +81-75-662-8967 FAX +81-75-662-8973

